

## Mouse Arrb2/Beta-arrestin-2 ELISA Kit

Catalog No: #EK19887

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Arrb2/Beta-arrestin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-arrestin-2
Target Name	Beta-arrestin-2
Other Names	Arrestin beta-2
Accession No.	Q91YI4
Uniprot	Q91YI4
GeneID	216869;

## Application Details

Assay Type: Sandwich

Detect Range:0.156-10 ng/mL Sensitivity:less than 0.081 ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm

gene\_full\_name:Arrb2

Note: This product is for in vitro research use only