

Rat Arrb2/Beta-arrestin-2 ELISA Kit

Catalog No: #EK19888

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Rat Arrb2/Beta-arrestin-2 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Beta-arrestin-2 |
| Target Name | Beta-arrestin-2 |
| Other Names | Arrestin beta-2 |
| Accession No. | P29067 |
| Uniprot | P29067 |
| GeneID | 25388; |

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Arrb2

Note: This product is for in vitro research use only