

Mouse Apoc2/Apolipoprotein C-II ELISA Kit

Catalog No: #EK19889

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apoc2/Apolipoprotein C-II ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apolipoprotein C-II
Target Name	Apolipoprotein C-II
Other Names	Apo-CII, Apolipoprotein C2
Accession No.	Q05020
Uniprot	Q05020
GeneID	105886299;11813;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.186ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apoc2

Note: This product is for in vitro research use only