

Pig CFB/Complement factor B ELISA Kit

Catalog No: #EK19891



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pig CFB/Complement factor B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Complement factor B
Target Name	Complement factor B
Other Names	C3/C5 convertase, Properdin factor B, BF
Accession No.	Q03710
Uniprot	Q03710

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.049ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CFB

Note: This product is for in vitro research use only