

Rat Aldh2/Aldehyde dehydrogenase, mitochondrial ELISA Kit

Catalog No: #EK19892

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Aldh2/Aldehyde dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aldehyde dehydrogenase, mitochondrial
Target Name	Aldehyde dehydrogenase, mitochondrial
Other Names	ALDH class 2,ALDH-E2,ALDH1
Accession No.	P11884
Uniprot	P11884
GeneID	29539;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.058ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Aldh2

Note: This product is for in vitro research use only