

# Mouse Cd80/T-lymphocyte activation antigen CD80 ELISA Kit

Catalog No: #EK19893

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cd80/T-lymphocyte activation antigen CD80 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-lymphocyte activation antigen CD80
Target Name	T-lymphocyte activation antigen CD80
Other Names	Activation B7-1 antigen,B7,CD80,B7
Accession No.	Q00609
Uniprot	Q00609
GeneID	12519;

## Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.189ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Cd80

Note: This product is for in vitro research use only