

## Chicken CKB/Creatine kinase B-type ELISA Kit

Catalog No: #EK19894

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken CKB/Creatine kinase B-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Creatine kinase B-type
Target Name	Creatine kinase B-type
Other Names	B-CK,Creatine kinase B chain
Accession No.	P05122
Uniprot	P05122
GeneID	396248;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CKB

Note: This product is for in vitro research use only