

## Dog CKB/Creatine kinase B-type ELISA Kit

Catalog No: #EK19895

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CKB/Creatine kinase B-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Creatine kinase B-type
Target Name	Creatine kinase B-type
Other Names	B-CK,Creatine kinase B chain
Accession No.	P05124
Uniprot	P05124
GeneID	100855552;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CKB

Note: This product is for in vitro research use only