

Mouse Ckb/Creatine kinase B-type ELISA Kit

Catalog No: #EK19896

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ckb/Creatine kinase B-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Creatine kinase B-type
Target Name	Creatine kinase B-type
Other Names	B-CK,Creatine kinase B chain,Ckbb
Accession No.	Q04447
Uniprot	Q04447
GeneID	12709;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ckb

Note: This product is for in vitro research use only