

Human BIRC5/Baculoviral IAP repeat-containing protein 5 ELISA Kit



Catalog No: #EK19897

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BIRC5/Baculoviral IAP repeat-containing protein 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Baculoviral IAP repeat-containing protein 5
Target Name	Baculoviral IAP repeat-containing protein 5
Other Names	Apoptosis inhibitor 4, Apoptosis inhibitor survivin, API4, IAP4
Accession No.	O15392
Uniprot	O15392
GeneID	332;

Application Details

Assay Type: Sandwich

Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul

Assay Time: 3-5h Detection wavelength: 450 nm

gene_full_name: BIRC5

Note: This product is for in vitro research use only