

Rat Birc5/Baculoviral IAP repeat-containing protein 5 ELISA Kit

Catalog No: #EK19899

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Birc5/Baculoviral IAP repeat-containing protein 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Baculoviral IAP repeat-containing protein 5
Target Name	Baculoviral IAP repeat-containing protein 5
Other Names	Apoptosis inhibitor survivin
Accession No.	Q9JHY7
Uniprot	Q9JHY7
GeneID	64041;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.64pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Birc5

Note: This product is for in vitro research use only