

Rat Cma1/Chymase ELISA Kit

Catalog No: #EK19901

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cma1/Chymase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chymase
Target Name	Chymase
Other Names	Alpha-chymase,Mast cell protease 3,rMCP-3,Mast cell protease 5,rMCP-5,Mast cell protease III,rMCP-III,Mcpt3
Accession No.	P50339
Uniprot	P50339
GeneID	25627;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.177ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cma1

Note: This product is for in vitro research use only