

Rat Artn/Artemin ELISA Kit

Catalog No: #EK19902



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                    |                                     |
|--------------------|-------------------------------------|
| Product Name       | Rat Artn/Artemin ELISA Kit          |
| Brief Description  | ELISA Kits                          |
| Applications       | ELISA                               |
| Species Reactivity | Rt                                  |
| Specificity        | Natural and recombinant Rat Artemin |
| Target Name        | Artemin                             |
| Accession No.      | Q6AYE8                              |
| Uniprot            | Q6AYE8                              |
| GeneID             | 362572;                             |

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Artn

Note: This product is for in vitro research use only