

Human C4BPB/C4b-binding protein beta chain ELISA Kit



Catalog No: #EK19903

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human C4BPB/C4b-binding protein beta chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C4b-binding protein beta chain
Target Name	C4b-binding protein beta chain
Accession No.	P20851
Uniprot	P20851
GeneID	725;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.18ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C4BPB

Note: This product is for in vitro research use only