

Human CCR3/C-C chemokine receptor type 3 ELISA Kit



Catalog No: #EK19905

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCR3/C-C chemokine receptor type 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C chemokine receptor type 3
Target Name	C-C chemokine receptor type 3
Other Names	C-C CKR-3, Eosinophil eotaxin receptor, CD193, CMKBR3
Accession No.	P51677
Uniprot	P51677
GeneID	1232;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCR3

Note: This product is for in vitro research use only