

Mouse Ccr3/Probable C-C chemokine receptor type 3 ELISA Kit

Catalog No: #EK19906

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccr3/Probable C-C chemokine receptor type 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Probable C-C chemokine receptor type 3
Target Name	Probable C-C chemokine receptor type 3
Other Names	C-C CKR-3, Macrophage inflammatory protein 1-alpha receptor-like 2, MIP-1 alpha RL2, CD193, Cmkbr1l2, Cmkbr3
Accession No.	P51678
Uniprot	P51678
GeneID	12771;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccr3

Note: This product is for in vitro research use only