

Mouse Adamts2/A disintegrin and metalloproteinase with thrombospondin motifs 2 ELISA Kit

Catalog No: #EK19908

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adamts2/A disintegrin and metalloproteinase with thrombospondin motifs 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 2
Target Name	A disintegrin and metalloproteinase with thrombospondin motifs 2
Other Names	ADAM-TS 2,Procollagen I N-proteinase,PC I-NP,Procollagen I/II amino propeptide-processing enzyme,Procollagen N-endopeptidase,pNPI
Accession No.	Q8C9W3
Uniprot	Q8C9W3
GeneID	216725;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.412ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Adamts2

Note: This product is for in vitro research use only