

Human AKR1C2/Aldo-keto reductase family 1 member C2 ELISA Kit

Catalog No: #EK19909

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AKR1C2/Aldo-keto reductase family 1 member C2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aldo-keto reductase family 1 member C2
Target Name	Aldo-keto reductase family 1 member C2
Other Names	3-alpha-HSD3, Chlordecone reductase homolog HAKRD, Dihydrodiol dehydrogenase 2, DD-2, Dihydrodiol dehydrogenase/bile acid-binding protein, DD/BABP, Trans-1,2-dihydrobenzene-1,2-diol dehydrogenase, Type III 3-alpha-hydroxysteroid dehydrogenase, DDH2
Accession No.	P52895
Uniprot	P52895
GeneID	1646;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.15ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
 Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AKR1C2

Note: This product is for in vitro research use only