

Bovine C1QTNF9/Complement C1q and tumor necrosis factor-related protein 9 ELISA Kit



Catalog No: #EK19910

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine C1QTNF9/Complement C1q and tumor necrosis factor-related protein 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Complement C1q and tumor necrosis factor-related protein 9
Target Name	Complement C1q and tumor necrosis factor-related protein 9
Accession No.	Q0II24
Uniprot	Q0II24
GeneID	523616;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C1QTNF9

Note: This product is for in vitro research use only