

Human AIM2/Interferon-inducible protein AIM2 ELISA Kit



Catalog No: #EK19912

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AIM2/Interferon-inducible protein AIM2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Interferon-inducible protein AIM2
Target Name	Interferon-inducible protein AIM2
Other Names	Absent in melanoma 2
Accession No.	O14862
Uniprot	O14862
GeneID	9447;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 24pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AIM2

Note: This product is for in vitro research use only