

Rat Csn1s1/Alpha-S1-casein ELISA Kit

Catalog No: #EK19914

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Csn1s1/Alpha-S1-casein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-S1-casein
Target Name	Alpha-S1-casein
Other Names	Alpha-casein,Csna
Accession No.	P02661
Uniprot	P02661
GeneID	24284;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Csn1s1

Note: This product is for in vitro research use only