

# Mouse Crhbp/Corticotropin-releasing factor-binding protein ELISA Kit

Catalog No: #EK19916

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Crhbp/Corticotropin-releasing factor-binding protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Corticotropin-releasing factor-binding protein
Target Name	Corticotropin-releasing factor-binding protein
Other Names	CRF-BP, Corticotropin-releasing hormone-binding protein, CRH-BP
Accession No.	Q60571
Uniprot	Q60571
GeneID	12919;

## Application Details

Assay Type: Sandwich Detect Range: 6.25-400ng/mL Sensitivity: 1.56ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Crhbp

Note: This product is for in vitro research use only