

## Mouse Apom/Apolipoprotein M ELISA Kit

Catalog No: #EK19918

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Apom/Apolipoprotein M ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apolipoprotein M
Target Name	Apolipoprotein M
Other Names	Apo-M, Ng20
Accession No.	Q9Z1R3
Uniprot	Q9Z1R3
GeneID	55938;

## Application Details

Assay Type: Sandwich Detect Range: 7.8-500ng/mL Sensitivity: 4.09ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Apom

Note: This product is for in vitro research use only