

Human APOD/Apolipoprotein D ELISA Kit

Catalog No: #EK19919

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human APOD/Apolipoprotein D ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Apolipoprotein D
Target Name	Apolipoprotein D
Other Names	Apo-D
Accession No.	P05090
Uniprot	P05090
GeneID	347;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.089 ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: APOD

Note: This product is for in vitro research use only