

Mouse Apod/Apolipoprotein D ELISA Kit

Catalog No: #EK19920



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Apod/Apolipoprotein D ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apolipoprotein D
Target Name	Apolipoprotein D
Other Names	Apo-D
Accession No.	P51910
Uniprot	P51910
GeneID	11815;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.092ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Apod

Note: This product is for in vitro research use only