

Rat Apod/Apolipoprotein D ELISA Kit

Catalog No: #EK19921

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Apod/Apolipoprotein D ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Apolipoprotein D
Target Name	Apolipoprotein D
Other Names	Apo-D
Accession No.	P23593
Uniprot	P23593
GeneID	25239;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.16ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Apod

Note: This product is for in vitro research use only