

Mouse Amy2/Pancreatic alpha-amylase ELISA Kit

Catalog No: #EK19922

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Amy2/Pancreatic alpha-amylase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Pancreatic alpha-amylase
Target Name	Pancreatic alpha-amylase
Other Names	PA,1,4-alpha-D-glucan glucanohydrolase,Amy2a
Accession No.	P00688
Uniprot	P00688
GeneID	100043684;100043686;100043688;109959;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.164ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Amy2

Note: This product is for in vitro research use only