

## Rat Amy2/Pancreatic alpha-amylase ELISA Kit

Catalog No: #EK19923

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Amy2/Pancreatic alpha-amylase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Pancreatic alpha-amylase
Target Name	Pancreatic alpha-amylase
Other Names	PA,1,4-alpha-D-glucan glucanohydrolase
Accession No.	P00689
Uniprot	P00689
GeneID	497039;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Amy2

Note: This product is for in vitro research use only