

Mouse Bst1/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2 ELISA Kit

Catalog No: #EK19924

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bst1/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2
Target Name	ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2
Other Names	ADP-ribosyl cyclase 2, Antigen BP3, BP-3 alloantigen, Bone marrow stromal antigen 1, BST-1, Cyclic ADP-ribose hydrolase 2, cADPr hydrolase 2, Leukocyte antigen 65, Ly-65, CD157, Bp-3, Bp3, Ly65
Accession No.	Q64277
Uniprot	Q64277
GeneID	12182;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.088ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bst1

Note: This product is for in vitro research use only