

Rat Bst1/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2 ELISA Kit

Catalog No: #EK19925

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bst1/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2
Target Name	ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2
Other Names	ADP-ribosyl cyclase 2, Bone marrow stromal antigen 1, BST-1, Cyclic ADP-ribose hydrolase 2, cADPr hydrolase 2, CD157
Accession No.	Q63072
Uniprot	Q63072
GeneID	81506;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bst1

Note: This product is for in vitro research use only