

Human AIMP1/Aminoacyl tRNA synthase complex-interacting multifunctional protein 1 ELISA Kit



Catalog No: #EK19928

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AIMP1/Aminoacyl tRNA synthase complex-interacting multifunctional protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aminoacyl tRNA synthase complex-interacting multifunctional protein 1
Target Name	Aminoacyl tRNA synthase complex-interacting multifunctional protein 1
Other Names	Multisynthase complex auxiliary component p43,EMAP2,SCYE1
Accession No.	Q12904
Uniprot	Q12904
GeneID	9255;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/ml Sensitivity: 0.421ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AIMP1

Note: This product is for in vitro research use only