

Human CSF3R/Granulocyte colony-stimulating factor receptor ELISA Kit

Catalog No: #EK19929

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CSF3R/Granulocyte colony-stimulating factor receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Granulocyte colony-stimulating factor receptor
Target Name	Granulocyte colony-stimulating factor receptor
Other Names	G-CSF receptor,CD114,GCSFR
Accession No.	Q99062
Uniprot	Q99062
GeneID	1441;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.39ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CSF3R

Note: This product is for in vitro research use only