

Rat Aqp9/Aquaporin-9 ELISA Kit

Catalog No: #EK19931

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Aqp9/Aquaporin-9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aquaporin-9
Target Name	Aquaporin-9
Other Names	AQP-9, Aquaglyceroporin-9
Accession No.	P56627
Uniprot	P56627
GeneID	65054;

Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.68ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Aqp9

Note: This product is for in vitro research use only