

Mouse Adipor1/Adiponectin receptor protein 1 ELISA Kit



Catalog No: #EK19932

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adipor1/Adiponectin receptor protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adiponectin receptor protein 1
Target Name	Adiponectin receptor protein 1
Accession No.	Q91VH1
Uniprot	Q91VH1
GeneID	72674;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Adipor1

Note: This product is for in vitro research use only