

CACNG8 Antibody

Catalog No: #44701

Package Size: #44701-1 50ul #44701-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CACNG8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CACNG8 antibody detects endogenous levels of total CACNG8
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CACNG8
Target Name	CACNG8
Other Names	CACNG antibody CACNG6 antibody Cacng8 antibody Calcium channel voltage dependent gamma subunit 8 antibody CCG8_HUMAN antibody Neuronal voltage gated calcium channel gamma 8 subunit antibody Neuronal voltage-gated calcium channel gamma-8 subunit antibody TARP gamma-8 antibody Transmembrane AMPAR regu
Accession No.	Swiss-Prot#: Q8WXS5 NCBI Gene ID:59283
Uniprot	Q8WXS5
GeneID	59283;
Calculated MW	43 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only