

Mouse Bmx/Cytoplasmic tyrosine-protein kinase BMX ELISA Kit



Catalog No: #EK19933

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bmx/Cytoplasmic tyrosine-protein kinase BMX ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cytoplasmic tyrosine-protein kinase BMX
Target Name	Cytoplasmic tyrosine-protein kinase BMX
Other Names	Bone marrow tyrosine kinase gene in chromosome X protein homolog
Accession No.	P97504
Uniprot	P97504
GeneID	12169;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bmx

Note: This product is for in vitro research use only