

Rat Azgp1/Zinc-alpha-2-glycoprotein ELISA Kit

Catalog No: #EK19934

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Azgp1/Zinc-alpha-2-glycoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Zinc-alpha-2-glycoprotein
Target Name	Zinc-alpha-2-glycoprotein
Other Names	Zn-alpha-2-GP
Accession No.	Q63678
Uniprot	Q63678

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.164ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Azgp1

Note: This product is for in vitro research use only