

## Mouse Capn1/Calpain-1 catalytic subunit ELISA Kit

Catalog No: #EK19935

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Capn1/Calpain-1 catalytic subunit ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calpain-1 catalytic subunit
Target Name	Calpain-1 catalytic subunit
Other Names	Calcium-activated neutral proteinase 1,CANP 1,Calpain mu-type,Calpain-1 large subunit,Micromolar-calpain,muCANP,Canp1,Capa1
Accession No.	O35350
Uniprot	O35350
GeneID	12333;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.168ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Capn1

Note: This product is for in vitro research use only