

Human CAPN3/Calpain-3 ELISA Kit

Catalog No: #EK19937

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CAPN3/Calpain-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calpain-3
Target Name	Calpain-3
Other Names	Calcium-activated neutral proteinase 3,CANP 3,Calpain L3,Calpain p94,Muscle-specific calcium-activated neutral protease 3,New calpain 1,nCL-1,CANP3,CANPL3,NCL1
Accession No.	P20807
Uniprot	P20807
GeneID	825;

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.397ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CAPN3

Note: This product is for in vitro research use only