

Mouse Capn3/Calpain-3 ELISA Kit

Catalog No: #EK19938

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Capn3/Calpain-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calpain-3
Target Name	Calpain-3
Other Names	Calcium-activated neutral proteinase 3,CANP 3,Calpain L3,Calpain p94,Muscle-specific calcium-activated neutral protease 3
Accession No.	Q64691
Uniprot	Q64691
GeneID	12335;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.79ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Capn3

Note: This product is for in vitro research use only