

## Human CAPN5/Calpain-5 ELISA Kit

Catalog No: #EK19939

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CAPN5/Calpain-5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calpain-5
Target Name	Calpain-5
Other Names	Calpain htra-3, New calpain 3, nCL-3, NCL3
Accession No.	O15484
Uniprot	O15484
GeneID	726;

## Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.791ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CAPN5

Note: This product is for in vitro research use only