

Mouse Capn5/Calpain-5 ELISA Kit

Catalog No: #EK19940

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Capn5/Calpain-5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calpain-5
Target Name	Calpain-5
Other Names	New calpain 3,nCL-3,Ncl3
Accession No.	O08688
Uniprot	O08688
GeneID	12337;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.792ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Capn5

Note: This product is for in vitro research use only