

## Mouse Capn6/Calpain-6 ELISA Kit

Catalog No: #EK19941

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Capn6/Calpain-6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calpain-6
Target Name	Calpain-6
Other Names	Capa6
Accession No.	O35646
Uniprot	O35646
GeneID	12338;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Capn6

Note: This product is for in vitro research use only