

## Human ACY1/Aminoacylase-1 ELISA Kit

Catalog No: #EK19942

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ACY1/Aminoacylase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aminoacylase-1
Target Name	Aminoacylase-1
Other Names	ACY-1, N-acyl-L-amino-acid amidohydrolase
Accession No.	Q03154
Uniprot	Q03154
GeneID	95;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.399ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ACY1

Note: This product is for in vitro research use only