

## Mouse Acy1/Aminoacylase-1 ELISA Kit

Catalog No: #EK19943

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Acy1/Aminoacylase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Aminoacylase-1
Target Name	Aminoacylase-1
Other Names	ACY-1,N-acyl-L-amino-acid amidohydrolase
Accession No.	Q99JW2
Uniprot	Q99JW2
GeneID	109652;

## Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.325ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Acy1

Note: This product is for in vitro research use only