

Pig ACY1/Aminoacylase-1 ELISA Kit

Catalog No: #EK19944

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig ACY1/Aminoacylase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Aminoacylase-1
Target Name	Aminoacylase-1
Other Names	ACY-1, N-acyl-L-amino-acid amidohydrolase
Accession No.	P37111
Uniprot	P37111
GeneID	396930;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.176ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ACY1

Note: This product is for in vitro research use only