

Rat Acy1a/Aminoacylase-1A ELISA Kit

Catalog No: #EK19945

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Acy1a/Aminoacylase-1A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aminoacylase-1A
Target Name	Aminoacylase-1A
Other Names	ACY-1A,ACY IA,N-acyl-L-amino-acid amidohydrolase,Acy1
Accession No.	Q6AYS7
Uniprot	Q6AYS7
GeneID	300981;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Acy1a

Note: This product is for in vitro research use only