

Mouse Aspa/Aspartoacylase ELISA Kit

Catalog No: #EK19946

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aspa/Aspartoacylase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Aspartoacylase
Target Name	Aspartoacylase
Other Names	Aminoacylase-2,ACY-2
Accession No.	Q8R3P0
Uniprot	Q8R3P0
GeneID	11484;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aspa

Note: This product is for in vitro research use only