

## Rat Aspa/Aspartoacylase ELISA Kit

Catalog No: #EK19947

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Aspa/Aspartoacylase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aspartoacylase
Target Name	Aspartoacylase
Other Names	Aminoacylase-2,ACY-2
Accession No.	Q9R1T5
Uniprot	Q9R1T5
GeneID	79251;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.091ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Aspa

Note: This product is for in vitro research use only