

## Mouse Casp7/Caspase-7 ELISA Kit

Catalog No: #EK19948

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Casp7/Caspase-7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Caspase-7
Target Name	Caspase-7
Other Names	CASP-7, Apoptotic protease Mch-3, Cysteine protease LICE2, Lice2, Mch3
Accession No.	P97864
Uniprot	P97864
GeneID	12369;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.092ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Casp7

Note: This product is for in vitro research use only